

QUALITY CONTROL OF PROTEOMICS DATA FROM LIQUID CHROMATOGRAPHY - MASS SPECTROMETER

Submitted by

Kishore Tumarada

Table of Contents

| | | |
|----|--|----|
| 1. | ABSTRACT..... | 3 |
| 2. | INTRODUCTION (BACKGROUND OF THE INTERNSHIP PROJECT, AND THE MOTIVATION OF THE STUDY) | 3 |
| 3. | DATA AND METHODS..... | 5 |
| 4. | RESULTS AND DISCUSSION | 9 |
| 5. | CONCLUSION | 21 |
| 6. | REFERENCES | 21 |
| 7. | TABLE AND FIGURES OF DATA ANALYSIS | 22 |
| 8. | APPENDIX – PYTHON CODE | 24 |

1. Abstract

- 1.1. Quality control of Liquid Chromatography - Mass Spectrometer system is essential for appropriate statistical analysis of output data. To monitor the system performance, a set of standard peptides and fragment ions have been extracted based on conditions of maximum intensity, uniformly sampled across whole retention time, and lowest retention time CV. In particular, a total of 21 peptides have been extracted along with 5 fragment ions with non-neutral losses. These standards are used as the basis for analyzing successive LC-MS outputs for different Raw files over a period of time. Finally, visualizations of various parameters - m/z , intensity, resolution, reduced mobility, ccs, retention time, retention length (FWHM), ion mobility index length (FWHM), fragment ion intensities - have been created to see the pattern of peptides and monitor the performance of the LC-MS system.

2. Introduction (background of the internship project, and the motivation of the study)

- 2.1. Proteomics refers to the large-scale experimental analysis of proteins and proteomes. It has enabled identification of ever-increasing number of proteins. It covers exploration of proteomes from overall level of protein composition, structure, and activity.
- 2.2. It is a crucial domain in modern biological and biomedical research. Presently, Liquid chromatography (LC) followed by Mass spectrometry (MS) is the preferred method to identify and quantify complex protein samples. The importance of these techniques is demonstrated by their use in large-scale research initiatives, such as the ongoing Human Proteome Project (HPP).




Figure 1 Workflow of LC-MS system

2.3. Figure 1 shows the workflow of LC-MS system. Firstly, proteins in a biological sample are digested into peptides by tryptic digestion. In the next stage, peptides are separated based on retention time and ion mobility. In the last stage, based on spectrum match scoring and database search, proteins are identified.




Figure 2. A diagram highlighting some of the major sources of variability in each of the stages of an LC-MS experiment.

2.4. However, the results of an experiment can still be subject to significant variability. This variability can originate from multiple sources as shown in Figure 2. Hence quality control is key for proteomic workflow.

2.5. Goals of this project are three-fold. Firstly, extract a series of peptides to serve as standards. Secondly, extract information only for this group of peptides from QC (Quality control) runs obtained periodically. Lastly, visualize the parameters of each of these standard peptides for subsequent QC runs overtime to monitor the performance of this system.

3. Data and methods

3.1. LC-MS output Raw data is high-resolution data without data labeling. MaxQuant, a quantitative proteomics software package, is used for analyzing and data labeling the raw data. The output is used as input for this project.

3.2. The following datasets from MaxQuant output are considered for creating benchmark peptides:

3.2.1. **Allpeptides file** - all detected LC-MS features (feature characteristics; can be e.g. used to plot all features against the ones targeted for MS/MS and identified ones)

| | Raw file | Charge | m/z | Mass | Resolution | Number of data points | Number of frames | Number of isotopic peaks | Isotope correlation | Mass fractional part | ... |
|---|------------------------------------|--------|-----------|-----------|--------------|-----------------------|------------------|--------------------------|---------------------|----------------------|-----|
| 0 | HeLa-iRT-200ng-90min_Slot1-3_01_63 | 1 | 368.42412 | 367.41684 | 23132.500394 | 91 | 5 | 2 | 0.998080 | 0.416839 | ... |
| 1 | HeLa-iRT-200ng-90min_Slot1-3_01_63 | 1 | 367.18877 | 366.18150 | 22788.935224 | 2394 | 20 | 2 | 0.998260 | 0.181498 | ... |
| 2 | HeLa-iRT-200ng-90min_Slot1-3_01_63 | 1 | 369.18406 | 368.17679 | 20780.897048 | 726 | 13 | 2 | 0.993305 | 0.176788 | ... |
| 3 | HeLa-iRT-200ng-90min_Slot1-3_01_63 | 1 | 339.15809 | 338.15082 | 23096.147644 | 389 | 10 | 2 | 0.993477 | 0.150818 | ... |
| 4 | HeLa-iRT-200ng-90min_Slot1-3_01_63 | 1 | 351.37303 | 350.36576 | 22956.937759 | 258 | 7 | 2 | 0.999412 | 0.365758 | ... |

5 rows × 24 columns

| ... | Min frame index | Max frame index | Ion mobility index | Ion mobility index length | Ion mobility index length (FWHM) | Intensity | Intensities | Number of pasef MS/MS | Pasef MS/MS IDs | MS/MS scan number |
|-----|-----------------|-----------------|--------------------|---------------------------|----------------------------------|-----------|-------------|-----------------------|-----------------|-------------------|
| ... | 4965 | 4969 | 654 | 18 | 15 | 1498.6 | NaN | 0 | NaN | NaN |
| ... | 5088 | 5107 | 669 | 33 | 24 | 17524.0 | NaN | 0 | NaN | NaN |
| ... | 5089 | 5101 | 669 | 27 | 18 | 7609.5 | NaN | 0 | NaN | NaN |
| ... | 4959 | 4968 | 714 | 30 | 21 | 8555.6 | NaN | 0 | NaN | NaN |
| ... | 4765 | 4771 | 717 | 24 | 18 | 6640.4 | NaN | 0 | NaN | NaN |

Table 1. Top 10 rows of Allpeptides dataset

3.2.2. Msms file: identified MS/MS events

| | Raw file | Scan number | Scan index | Sequence | Length | Missed cleavages | Modifications | Modified sequence | Oxidation (M) Probabilities | Oxidation (M) Score diffs | ... |
|---|------------------------------------|-------------|------------|------------------|--------|------------------|-------------------------|---|-----------------------------|---------------------------|-----|
| 0 | HeLa-iRT-200ng-90min_Slot1-3_01_63 | 26150 | 26150 | AAAAAAAAAAGAAGGR | 16 | 0 | Acetyl (Protein N-term) | _(Acetyl (Protein N-term))AAAAAAAAAAGAAGGR_ | NaN | NaN | ... |
| 1 | HeLa-iRT-200ng-90min_Slot1-3_01_64 | 24634 | 24634 | AAAAAAAAAAGAAGGR | 16 | 0 | Acetyl (Protein N-term) | _(Acetyl (Protein N-term))AAAAAAAAAAGAAGGR_ | NaN | NaN | ... |
| 2 | HeLa-iRT-200ng-90min_Slot1-3_01_65 | 24700 | 24700 | AAAAAAAAAAGAAGGR | 16 | 0 | Acetyl (Protein N-term) | _(Acetyl (Protein N-term))AAAAAAAAAAGAAGGR_ | NaN | NaN | ... |
| 3 | HeLa-iRT-100ng-90min_Slot1-3_01_61 | 17560 | 17560 | AAAAAAAAAAGAAGGR | 16 | 0 | Acetyl (Protein N-term) | _(Acetyl (Protein N-term))AAAAAAAAAAGAAGGR_ | NaN | NaN | ... |
| 4 | HeLa-iRT-100ng-90min_Slot1-3_01_62 | 18039 | 18039 | AAAAAAAAAAGAAGGR | 16 | 0 | Acetyl (Protein N-term) | _(Acetyl (Protein N-term))AAAAAAAAAAGAAGGR_ | NaN | NaN | ... |

5 rows × 63 columns

| ... | All sequences | All modified sequences | Reporter PIF | Reporter fraction | id | Protein group IDs | Peptide ID | Mod. peptide ID | Evidence ID | Oxidation (M) site IDs |
|-----|--|--|--------------|-------------------|----|-------------------|------------|-----------------|-------------|------------------------|
| ... | AAAAAAAAAAGAAGGR;AAAAAETPEVLR;PNLSGIPGESNR | _(Acetyl (Protein N-term))AAAAAAAAAAGAAGGR_; | NaN | NaN | 0 | 2479 | 0 | 0 | 0 | NaN |
| ... | AAAAAAAAAAGAAGGR;AAAAAETPEVLR;QPSRQSERPR | _(Acetyl (Protein N-term))AAAAAAAAAAGAAGGR_; | NaN | NaN | 1 | 2479 | 0 | 0 | 1 | NaN |
| ... | AAAAAAAAAAGAAGGR;AAAAAETPEVLR;IQRATQEPVAK | _(Acetyl (Protein N-term))AAAAAAAAAAGAAGGR_; | NaN | NaN | 2 | 2479 | 0 | 0 | 2 | NaN |
| ... | AAAAAAAAAAGAAGGR;AAAAAETPEVLR;IVMNRNNVHK | _(Acetyl (Protein N-term))AAAAAAAAAAGAAGGR_; | NaN | NaN | 3 | 2479 | 0 | 0 | 3 | NaN |
| ... | AAAAAAAAAAGAAGGR;AAAAAETPEVLR;RKNQSGTMFR | _(Acetyl (Protein N-term))AAAAAAAAAAGAAGGR_; | NaN | NaN | 4 | 2479 | 0 | 0 | 4 | NaN |

Table 2. Top 10 rows of msms dataset

3.2.3. **Evidence file:** all identified LC-MS-features, identified by MS/MS or by matching between runs, MS/MS events per identified feature, scores, mass deviation, feature information, useful table also for troubleshooting.

| | Sequence | Length | Modifications | Modified sequence | Oxidation (M) Probabilities | Oxidation (M) Score Diffs | Acetyl (Protein N-term) | Oxidation (M) | Missed cleavages | Proteins | ... | R |
|---|-----------------|--------|-------------------------|--|-----------------------------|---------------------------|-------------------------|---------------|------------------|-----------------|-----|---|
| 0 | AAAAAAAAAGAAGGR | 16 | Acetyl (Protein N-term) | _(Acetyl (Protein N-term))AAAAAAAAAGAAGGR_ | NaN | NaN | 1 | 0 | 0 | Q86U42;Q86U42-2 | ... | |
| 1 | AAAAAAAAAGAAGGR | 16 | Acetyl (Protein N-term) | _(Acetyl (Protein N-term))AAAAAAAAAGAAGGR_ | NaN | NaN | 1 | 0 | 0 | Q86U42;Q86U42-2 | ... | |
| 2 | AAAAAAAAAGAAGGR | 16 | Acetyl (Protein N-term) | _(Acetyl (Protein N-term))AAAAAAAAAGAAGGR_ | NaN | NaN | 1 | 0 | 0 | Q86U42;Q86U42-2 | ... | |
| 3 | AAAAAAAAAGAAGGR | 16 | Acetyl (Protein N-term) | _(Acetyl (Protein N-term))AAAAAAAAAGAAGGR_ | NaN | NaN | 1 | 0 | 0 | Q86U42;Q86U42-2 | ... | |
| 4 | AAAAAAAAAGAAGGR | 16 | Acetyl (Protein N-term) | _(Acetyl (Protein N-term))AAAAAAAAAGAAGGR_ | NaN | NaN | 1 | 0 | 0 | Q86U42;Q86U42-2 | ... | |

5 rows x 71 columns

| ... | Reporter fraction | Reverse | Potential contaminant | id | Protein group IDs | Peptide ID | Mod. peptide ID | MS/MS IDs | Best MS/MS | Oxidation (M) site IDs |
|-----|-------------------|---------|-----------------------|----|-------------------|------------|-----------------|-----------|------------|------------------------|
| ... | NaN | NaN | NaN | 0 | 2479 | 0 | 0 | 0 | 0 | NaN |
| ... | NaN | NaN | NaN | 1 | 2479 | 0 | 0 | 1 | 1 | NaN |
| ... | NaN | NaN | NaN | 2 | 2479 | 0 | 0 | 2 | 2 | NaN |
| ... | NaN | NaN | NaN | 3 | 2479 | 0 | 0 | 3 | 3 | NaN |
| ... | NaN | NaN | NaN | 4 | 2479 | 0 | 0 | 4 | 4 | NaN |

Table 3 Top 10 rows of Evidence dataset

4. Results and discussion

4.1. Part 1 – Extraction of standard peptides

4.1.1. General principles for extracting Peptide sequences:

To pick these standard peptides, we need to generate a pool of suitable ones. Following are the general principles adopted to extract data:

4.1.1.1. Limit peptides with a charge of 2. Column “charge” in both msms.txt and evidence.txt.

4.1.1.2. Standard peptides should be uniformly distributed along the entire gradient, i.e., across total retention time taken by a solute to pass through the liquid chromatography column. For example, if the gradient is 120 minutes, we can choose peptide every 5 minutes, which corresponds to at least 24 peptides. The following figure shows the histogram for retention time across the entire data




Figure 3 Histogram of Retention time of all Raw files

4.1.1.3. As the same peptide can be found multiple times in multiple files, the retention time has a bit of variation. For this, we have evaluated the Coefficient of Variations (CV) and decide a proper CV value to filter the peptides. (The coefficient of variation (CV), which is calculated by dividing the standard deviation of peptide profiles by the mean, is reported as percentages.)

4.1.1.4. Peptides need to be present in every sample or Raw file. For instance, in the present data, MQ output is generated from 6 raw files. Table 4 shows the number of peptide sequences present in each raw file.

```
evidence_grouped.count()
```

[7]:

| Sequence | Raw file |
|------------------------------------|----------|
| HeLa-iRT-100ng-90min_Slot1-3_01_61 | 19069 |
| HeLa-iRT-100ng-90min_Slot1-3_01_62 | 19093 |
| HeLa-iRT-100ng-90min_Slot1-3_01_66 | 18443 |
| HeLa-iRT-200ng-90min_Slot1-3_01_63 | 22654 |
| HeLa-iRT-200ng-90min_Slot1-3_01_64 | 22641 |
| HeLa-iRT-200ng-90min_Slot1-3_01_65 | 22505 |

Table 4 number of peptide sequences across each raw file

- 4.1.1.5. Peptides need to have good signal intensity.
- 4.1.2. Based on these principles, out of all peptide sequences in every 5 minutes interval, peptides present in all raw files are filtered. From these peptides, those with the highest intensity and lowest Retention time CV is selected. As a result, a maximum of 24 peptides for 120 minutes retention time interval can be obtained.
- 4.1.3. However, in the selected MQ output, only 21 peptides met the conditions. Table 5 shows the distribution of these peptides across raw files:

| Raw file | HeLa-iRT-100ng-90min_Slot1-3_01_61 | HeLa-iRT-100ng-90min_Slot1-3_01_62 | HeLa-iRT-100ng-90min_Slot1-3_01_66 | HeLa-iRT-200ng-90min_Slot1-3_01_63 | HeLa-iRT-200ng-90min_Slot1-3_01_64 | HeLa-iRT-200ng-90min_Slot1-3_01_65 |
|--------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| Sequence | | | | | | |
| AAGVNVEPFWPGLFAK | 1 | 1 | 1 | 1 | 1 | 1 |
| APNTPASGANGDGMSQTQSGSTVK | 3 | 2 | 2 | 2 | 2 | 3 |
| EHALLAYTLGVK | 1 | 1 | 1 | 1 | 1 | 2 |
| EILVGDDVGQTVDDPYATFVK | 1 | 3 | 3 | 4 | 4 | 2 |
| ESTLHLVLR | 1 | 2 | 1 | 2 | 2 | 2 |
| FHVEEEGK | 1 | 1 | 1 | 1 | 1 | 1 |
| FMQISEDSTR | 2 | 2 | 1 | 1 | 1 | 3 |
| GTFIDPAAVIR | 4 | 5 | 2 | 3 | 4 | 2 |
| GYSFTTTER | 4 | 4 | 2 | 2 | 3 | 4 |
| IINEPTAAAIAYGLDK | 3 | 5 | 4 | 4 | 7 | 6 |
| ISVYYNEATGGK | 1 | 3 | 1 | 1 | 2 | 1 |
| LIAPVAEEETVPNNK | 1 | 4 | 1 | 2 | 3 | 3 |
| MSVQPTVSLGGFEITPPVLR | 2 | 2 | 2 | 2 | 2 | 2 |
| NHEEEMNALR | 2 | 1 | 1 | 1 | 1 | 3 |
| NTGICTIGPASR | 1 | 1 | 3 | 2 | 2 | 2 |
| QADTVYFLPITPQFVTEVIK | 1 | 2 | 2 | 1 | 1 | 1 |
| QGGLGPMNIPLVSDPK | 3 | 6 | 3 | 3 | 6 | 4 |
| TIGGGDDSFNTFFSETGAGK | 5 | 3 | 5 | 4 | 5 | 7 |
| VFLENVIR | 5 | 9 | 5 | 11 | 11 | 9 |
| VHGPQISGTTNKPKNK | 1 | 1 | 1 | 1 | 1 | 1 |
| YPIEHGIITNWDDMEK | 1 | 1 | 1 | 1 | 2 | 1 |

Table 5. Distribution of peptides across all raw files

4.1.4. **General principles for fragment ions extraction:**

These standard Peptides need to have good fragmentation ion signals (also called MS2 signals).

- 4.1.4.1. This info on fragmentation ions can be found in columns “Matches”, “Intensities”, and “Masses” in file msms file. Column “Matches” contains the types of fragmented ions, and column “Intensities” and “Masses” contain the intensity and mass of each ion in the column “Matches”, respectively.
- 4.1.4.2. Only list the type, intensity, and mass of the fragment ions above a certain intensity value, for example 200.
- 4.1.4.3. Also eliminate those with neutral losses, such as y_5-H_2O , or y_3-NH_3 , etc. only include those ions with types in the form of a letter (y or b) with a number (single or double digits).

4.1.5. Table 6 shows the top 5 rows of msms data, after filtering for standard peptides:

| id | Sequence | Raw file | Matches | Intensities | Masses |
|-----|------------------|------------------------------------|---|---|---|
| 742 | AAGVNVEPFWPGLFAK | HeLa-iRT-200ng-90min_Slot1-3_01_63 | y2;y3;y6;y7;y9;y10;y11;y12;y13;y14;y10-H2O;y13... | 20;1000;8000;2000;20000;9000;2000;400;200;20;1... | 218.155605126463;365.218184501103;632.37397515... |
| 743 | AAGVNVEPFWPGLFAK | HeLa-iRT-200ng-90min_Slot1-3_01_64 | y2;y3;y6;y7;y8;y9;y10;y11;y12;y13;y14;y10-H2O;... | 30;700;7000;2000;200;10000;6000;2000;2000;100;... | 218.152568074812;365.217331549252;632.37104926... |
| 744 | AAGVNVEPFWPGLFAK | HeLa-iRT-200ng-90min_Slot1-3_01_65 | y2;y3;y5;y6;y7;y8;y9;y7-NH3;y9-NH3;b3;b4;b5;b6... | 10;700;400;8000;2000;300;500;200;100;4;700;200... | 218.154799944179;365.217327494956;535.31939018... |
| 745 | AAGVNVEPFWPGLFAK | HeLa-iRT-100ng-90min_Slot1-3_01_61 | y3;y5;y6;y7;y8;y9;y10;y11;y12;y14;y10-H2O;y12-... | 90;50;1000;300;70;2000;1000;100;500;4;100;4;10... | 365.21733850014;535.318608615786;632.371543281... |
| 746 | AAGVNVEPFWPGLFAK | HeLa-iRT-100ng-90min_Slot1-3_01_62 | y3;y5;y6;y7;y8;y9;y10;y11;y12;y13;y14;y10-H2O;... | 300;200;3000;800;300;7000;3000;900;1000;20;5;4... | 365.218363026717;535.322899774525;632.37385499... |

Table 6. Top 5 rows of msms data after filtering standard peptides

4.1.6. From the Filtered msms data, Top 5 fragment ions with highest intensities, without non-neutral losses, for each of 21 peptide sequences have been extracted. Table 7 shows the top 10 rows of the final fragment ions table.

| | Sequence | matches | intensities | masses |
|-----|--------------------------|---------|-------------|-------------|
| 4 | AAGVNVEPFWPGLFAK | y9 | 20000 | 1062.573312 |
| 5 | AAGVNVEPFWPGLFAK | y10 | 9000 | 1191.617494 |
| 2 | AAGVNVEPFWPGLFAK | y6 | 8000 | 632.373975 |
| 12 | AAGVNVEPFWPGLFAK | b7 | 4000 | 641.322338 |
| 11 | AAGVNVEPFWPGLFAK | b6 | 4000 | 512.280951 |
| 39 | APNTPASGANGDGMSQTQSGSTVK | y21 | 1000 | 1966.860105 |
| 22 | APNTPASGANGDGMSQTQSGSTVK | b12 | 700 | 1053.454106 |
| 9 | APNTPASGANGDGMSQTQSGSTVK | y13 | 500 | 1297.600408 |
| 43 | APNTPASGANGDGMSQTQSGSTVK | b4 | 400 | 384.186522 |
| 160 | APNTPASGANGDGMSQTQSGSTVK | y23 | 400 | 2181.947540 |

Table 7. Top 10 rows of Processed msms data with 5 fragment ions for each standard peptide

4.1.7. Table 8 shows the final data for standard peptides, after extracting Retention length and Ion mobility index length (FWHM) from allpeptides file.

| | Sequence | mean_retention_time | cv_retention_time | mean_intensity | cv_intensity | mean_mz | mean_masserror_ppm |
|---|--------------------------|---------------------|-------------------|----------------|--------------|-------------|--------------------|
| 0 | AAGVNVEPFWPGLFAK | 94.475500 | 0.059395 | 540880.000000 | 65.552959 | 851.951215 | 0.451648 |
| 1 | AAGVNVEPFWPGLFAK | 94.475500 | 0.059395 | 540880.000000 | 65.552959 | 851.951215 | 0.451648 |
| 2 | AAGVNVEPFWPGLFAK | 94.475500 | 0.059395 | 540880.000000 | 65.552959 | 851.951215 | 0.451648 |
| 3 | AAGVNVEPFWPGLFAK | 94.475500 | 0.059395 | 540880.000000 | 65.552959 | 851.951215 | 0.451648 |
| 4 | AAGVNVEPFWPGLFAK | 94.475500 | 0.059395 | 540880.000000 | 65.552959 | 851.951215 | 0.451648 |
| 5 | APNTPASGANGDGMSQTQSGSTVK | 21.221071 | 17.391118 | 923627.428571 | 68.465606 | 1178.959899 | -0.188700 |
| 6 | APNTPASGANGDGMSQTQSGSTVK | 21.221071 | 17.391118 | 923627.428571 | 68.465606 | 1178.959899 | -0.188700 |
| 7 | APNTPASGANGDGMSQTQSGSTVK | 21.221071 | 17.391118 | 923627.428571 | 68.465606 | 1178.959899 | -0.188700 |
| 8 | APNTPASGANGDGMSQTQSGSTVK | 21.221071 | 17.391118 | 923627.428571 | 68.465606 | 1178.959899 | -0.188700 |
| 9 | APNTPASGANGDGMSQTQSGSTVK | 21.221071 | 17.391118 | 923627.428571 | 68.465606 | 1178.959899 | -0.188700 |

| mean_resolution | mean_1_k0 | mean_ccs | matches | intensities | masses | Retention length (FWHM) | Ion mobility index length (FWHM) |
|-----------------|-----------|------------|---------|-------------|-------------|-------------------------|----------------------------------|
| 22968.064642 | 1.005164 | 405.855600 | y9 | 20000 | 1062.573312 | 0.124333 | 36.000000 |
| 22968.064642 | 1.005164 | 405.855600 | y10 | 9000 | 1191.617494 | 0.124333 | 36.000000 |
| 22968.064642 | 1.005164 | 405.855600 | y6 | 8000 | 632.373975 | 0.124333 | 36.000000 |
| 22968.064642 | 1.005164 | 405.855600 | b7 | 4000 | 641.322338 | 0.124333 | 36.000000 |
| 22968.064642 | 1.005164 | 405.855600 | b6 | 4000 | 512.280951 | 0.124333 | 36.000000 |
| 23535.742750 | 1.121327 | 451.742764 | y21 | 1000 | 1966.860105 | 0.214357 | 43.071429 |
| 23535.742750 | 1.121327 | 451.742764 | b12 | 700 | 1053.454106 | 0.214357 | 43.071429 |
| 23535.742750 | 1.121327 | 451.742764 | y13 | 500 | 1297.600408 | 0.214357 | 43.071429 |
| 23535.742750 | 1.121327 | 451.742764 | b4 | 400 | 384.186522 | 0.214357 | 43.071429 |
| 23535.742750 | 1.121327 | 451.742764 | y23 | 400 | 2181.947540 | 0.214357 | 43.071429 |

Table 8. Top 10 rows of dataset with standard peptides and their fragment ions

4.1.8. Following are histograms of CV for Retention time and Intensities:




Figure 4 Histogram of CV for Retention time

Figure 4 shows that highest frequency of CV of Retention time (in %) is less than 3%.




Figure 5 Histogram of CV of intensity

Figure 5 shows that highest frequency of CV of intensity (in %) is just less than 100%, which indicates that there is very high variation in intensity for peptides.

4.2. **Part 2 – Processing of QC run files based on standard peptides**

4.2.1. In the second part of the project, the goal is to extract data related to each peptide and its corresponding 5 fragment ions for every date (encoded form), which is extracted from raw file name. Table 9 shows the top 10 rows related to standard peptides.

| | sequence | date_extracted | m/z | mass error [ppm] | intensity | resolution | 1/k0 | ccs | retention time | retention length (fwhm) | ion mobility index length (fwhm) |
|---|--------------------------|----------------|------------|------------------|------------|--------------|----------|------------|----------------|-------------------------|----------------------------------|
| 0 | AAGVNVEPFWPGLFAK | 2901 | 851.951215 | 1.627200 | 2598200.00 | 23450.621678 | 1.063992 | 429.608425 | 94.056000 | 0.117000 | 39.00 |
| 1 | AAGVNVEPFWPGLFAK | 2917 | 851.951215 | 2.324400 | 3687700.00 | 24007.566346 | 1.074078 | 433.681027 | 97.584000 | 0.077200 | 42.00 |
| 2 | AAGVNVEPFWPGLFAK | 2924 | 851.951215 | 0.739230 | 492210.00 | 23746.773100 | 1.063992 | 429.608424 | 97.630000 | 0.077200 | 39.00 |
| 3 | APNTPASGANGDGSMSTQSGSTVK | 2901 | 983.110501 | -0.221770 | 574315.00 | 22900.649950 | 1.093244 | 540.128871 | 18.865000 | 0.195000 | 39.00 |
| 4 | APNTPASGANGDGSMSTQSGSTVK | 2917 | 983.110501 | 0.444375 | 983192.50 | 22987.426177 | 1.091559 | 539.110465 | 26.571000 | 0.174000 | 39.00 |
| 5 | APNTPASGANGDGSMSTQSGSTVK | 2924 | 983.110501 | 1.673482 | 404172.75 | 22932.699314 | 1.096592 | 541.477693 | 26.621000 | 0.178750 | 49.50 |
| 6 | APNTPASGANGDGSMSTQSGSTVK | 2937 | 983.110501 | 1.511998 | 801507.50 | 22928.193979 | 1.087289 | 536.619982 | 27.238500 | 0.178750 | 39.75 |
| 7 | EHALLAYTLGVK | 2901 | 657.874445 | 2.083200 | 2729500.00 | 24488.932846 | 0.959590 | 388.376985 | 56.616000 | 0.253000 | 39.00 |
| 8 | EHALLAYTLGVK | 2917 | 584.889204 | -13.786267 | 2512370.00 | 24162.673736 | 0.891936 | 412.260072 | 64.845333 | 0.161333 | 28.00 |
| 9 | EHALLAYTLGVK | 2924 | 548.396584 | -1.822527 | 296831.50 | 24374.728500 | 0.861493 | 425.914839 | 64.743000 | 0.135000 | 24.00 |

Table 9. Top 10 rows of standard peptides info extracted cumulatively from successive QC runs data

4.2.2. Similarly, data related to fragment ions, with highest intensities, of each standard peptide for each date (encoded form) has also been extracted as shown in Table 10.

| | sequence | date_extracted | matches | intensities | masses |
|---|------------------|----------------|---------|-------------|-------------|
| 0 | AAGVNVEPFWPGLFAK | 2901 | y6 | 5000 | 632.376515 |
| 1 | AAGVNVEPFWPGLFAK | 2901 | y9 | 2000 | 1062.579134 |
| 2 | AAGVNVEPFWPGLFAK | 2901 | y10 | 200 | 1191.623330 |
| 3 | AAGVNVEPFWPGLFAK | 2901 | b6 | 2000 | 512.281990 |
| 4 | AAGVNVEPFWPGLFAK | 2901 | b7 | 2000 | 641.324449 |
| 5 | AAGVNVEPFWPGLFAK | 2917 | y6 | 4000 | 632.369478 |
| 6 | AAGVNVEPFWPGLFAK | 2917 | y9 | 1000 | 1062.564936 |
| 7 | AAGVNVEPFWPGLFAK | 2917 | y10 | 100 | 1191.609416 |
| 8 | AAGVNVEPFWPGLFAK | 2917 | b6 | 2000 | 512.277118 |
| 9 | AAGVNVEPFWPGLFAK | 2917 | b7 | 1000 | 641.318439 |

Table 10. Top 10 rows of fragment ions info of standard peptides extracted cumulatively from successive QC runs data

4.3. Part 3 – Visualization of QC runs data

- 4.3.1. In the last part, data is visualized to identify the changes in parameters over a period of time. For illustration, we have taken 3 peptide sequences to measure the change in parameters over a period of time.
- 4.3.2. Figure 6 shows that mass-over-charge ratio vs date for different peptides. It shows that m/z remains constant over a period of time for all 3 peptides.




Figure 6 Line plot of date vs mass-to-charge ratio

4.3.3. Figure 7 shows that ion mobility index length increasing on 3rd day and then falling for all 3 peptides.




Figure 7 Line plot of date vs ion mobility index length

4.3.4. Figure 8 shows that retention length decreases on the second day but again increases on next day. However, GYS-peptide (blue line) does not reach its previous length.




Figure 8 Line plot of date vs retention length

4.3.5. Figure 9 shows that retention time increases for all peptides. Ideally, if we use the same experimental conditions, the sample's retention time should not change significantly.




Figure 9 Line plot of date vs retention time

4.3.6. Figures 10 and 11 show that ion mobility ($1 / K_0$) and CCS (Collision Cross Section) of all three peptides have stable values.




Figure 10 Line plot of date vs inverse reduced mobility




Figure 11. Line plot of date vs ccs

4.3.7. Figure 12 shows that the resolution of APN-peptide is almost stable, whereas the resolution of GYS-peptide has very high variation over time.




Figure 12. Line plot of date vs resolution

4.3.8. Figures 13 and 14 also show that intensity and mass errors of GYS-peptide also have very high variation compared to the other two peptides.




Figure 13. Line plot of date vs intensity




Figure 14. Line plot of date vs mass error

4.3.9. Figures 15-17 show changes in intensities for the top 5 fragment ions (based on standard peptides data) for each of these peptides. For APN-peptide, only 4 out of 5 fragment ions are present in QC runs. Out of these intensities of b4 and b12 ions increase over time. In the case of GYS-peptide, intensities of all 5 ions fall on 2nd date and again raise after that. However, in the case of last NHE-peptide, ion intensities increase on the same day and decrease later.




Figure 15. Line plot of date vs fragment ions intensities for a selected peptide




Figure 16. Line plot of date vs fragment ions intensities for a selected peptide




Figure 17. Line plot of date vs fragment ions intensities for a selected peptide

- 4.4. In this way, we can look at the patterns of any of extracted standard peptides, and their corresponding fragment ions of all QC Runs data obtained by Mass spectrometer.
- 4.5. In summary, this whole workflow effectively monitors the performance of LC-MS system through the lens of selected standard peptides and aids in tuning the system for best performance.

5. Conclusion

- 5.1. LC-MS system is prone to many variabilities, which decrease the reliability of the system in identifying the proteins. To control the quality of the system, a series of peptides and fragment ions have been chosen as standards, and various parameters such as intensities, m/z , ion mobilities have been monitored for these standards across successive runs of the system. This would help in tuning different parameters of the LC-MS system to improve quality.

6. References

- 6.1. Quality control in mass spectrometry-based proteomics by Wout Bittremieux, David Tabb et al.
- 6.2. *pmartR*: Quality Control and Statistics for Mass Spectrometry-Based Biological Data Kelly G. Stratton, Bobbie-Jo M. Webb-Robertson et al.

7. Table and figures of data analysis

7.1. List of Figures:

| Title | Page |
|---|------|
| Figure 1. Workflow of LC-MS system | 4 |
| Figure 2. A diagram highlighting some of the major sources of variability in each of the stages of an LC-MS experiment. | 4 |
| Figure 3. Histogram of Retention time of all Raw files | 9 |
| Figure 4. Histogram of CV for Retention time | 14 |
| Figure 5. Histogram of CV of intensity | 14 |
| Figure 6. Line plot of date vs mass-to-charge ratio | 16 |
| Figure 7. Line plot of date vs ion mobility index length | 17 |
| Figure 8. Line plot of date vs retention length | 17 |
| Figure 9. Line plot of date vs retention time | 18 |
| Figure 10. Line plot of date vs inverse reduced mobility | 18 |
| Figure 11. Line plot of date vs ccs | 18 |
| Figure 12. Line plot of date vs resolution | 19 |
| Figure 13. Line plot of date vs intensity | 19 |
| Figure 14. Line plot of date vs mass error | 19 |
| Figure 15. Line plot of date vs fragment ions intensities for a selected peptide | 20 |
| Figure 16. Line plot of date vs fragment ions intensities for a selected peptide | 20 |
| Figure 17. Line plot of date vs fragment ions intensities for a selected peptide | 20 |

7.2. List of tables

| Title | Page |
|--|------|
| Table 1. Top 10 rows of Allpeptides dataset | 6 |
| Table 2. Top 10 rows of msms dataset | 7 |
| Table 3 Top 10 rows of Evidence dataset | 8 |
| Table 4 number of peptide sequences across each raw file | 10 |
| Table 5. Distribution of peptides across all raw files | 11 |
| Table 6. Top 5 rows of msms data after filtering standard peptides | 12 |
| Table 7. Top 10 rows of Processed msms data with 5 fragment ions for each standard peptide | 12 |
| Table 8. Top 10 rows of dataset with standard peptides and their fragment ions | 13 |
| Table 9. Top 10 rows of standard peptides info extracted cumulatively from successive QC runs data | 15 |
| Table 10. Top 10 rows of fragment ions info of standard peptides extracted cumulatively from successive QC runs data | 15 |

8. Appendix – Python code

8.1. Part 1

```
#!/usr/bin/env python
# coding: utf-8
# # Code to extract benchmark peptides from Spectrometer data
# In[1]:
import pandas as pd
import numpy as np
import os
from pandas import ExcelWriter
import seaborn as sns
import matplotlib.pyplot as plt
import re
import warnings
# # Reading all files in the test/txt directory
# # Please enter the location of Maxquant output in test_directory variable
# In[2]:
test_directory = r'/Volumes/Transcend/OneDrive/UoH DS/Summer 2020/Spectrometry
lab internship/new_MQoutput/MaxQuantOutput/test/txt'
# In[3]:
#display all files in the directory
for root, dirs, files in os.walk(test_directory, topdown=False):
    file_names = []
    print('root: ', root)
    print('filenames--')
    for name in files:
        if name[0] is not '.':
            print(os.path.join(root, name))
            file_names.append(os.path.join(root, name))

# In[4]:
warnings.filterwarnings("ignore")
#Reading filenames of evidence and msms files
for name in file_names:
    if name.endswith('evidence.txt'):
        file_evidence = name
    elif name.endswith('msms.txt'):
        file_msms = name
    elif name.endswith('allPeptides.txt'):
        file_allpeptides = name

# In[5]:
print(file_evidence)
```



```

df_evidence = pd.read_table(file_evidence)
df_evidence.head()

# In[6]:
print(file_msms)
df_msms = pd.read_table(file_msms)
df_msms.head()

# In[7]:
print(file_allpeptides)
df_allpeptides = pd.read_table(file_allpeptides)
df_allpeptides.head()

# ## Info needed for each peptide:
#     Peptide sequence - To be extracted from evidence file
#     m/z             - Evidence file
#     mass error (ppm) - Evidence file
#     Retention time (FWHM) - allpeptides
#     Retention length - Evidence file
#     Intensity       - Evidence file
#     Resolution      - Evidence file
#     1/k0            - Evidence file
#     CCS             - Evidence file
#     Ion mobility index length (FWHM) - allpeptides
#
# ## The following info needed for fragment ions:
#     Matches - msms file
#     Intensities- msms file
#     Masses- msms file

# # Evidence file processing

# ## Grouping content to Raw files

# In[8]:
evidence_grouped = df_evidence.groupby(['Raw file'])
evidence_grouped.count()

# In[9]:
evidence_grouped['Retention time'].agg(['mean', 'min', 'max'])

# ## Extracting evidence data for charge = 2
# In[10]:

```

```

#Filter for charge
df_evidence_q2 = df_evidence.query('Charge == 2')
df_evidence_q2.describe()

# ## Convert Retention time column to date time format
# In[11]:

df_evidence_q2['RT_timeform'] = pd.to_timedelta(df_evidence_q2['Retention time'],
unit='minutes')
df_evidence_q2.head()

# ## Resample for every 5 mins and extract peptides list
# In[12]:
#Time resampling of evidence files
resampled = df_evidence_q2.resample(rule = '5T',on = 'RT_timeform')
num_rawfiles = df_evidence_q2['Raw file'].unique().size

#aggregation function for grouped and resampled data
def fun_count_rawfiles(series):
    return series.unique().size

#aggregate dfunction for calculation of cv
def func_cv(df):
    return np.abs(df.std()*100/df.mean())

#creating empty list for saving standard peptides later
selected_peptides = []

#Read each group of 5 mins of peptides and extract 1 peptide from each
for name, group in resampled:
    print("Group: ", name)
    print("-" * 27)
    sample_group = group[['Sequence', 'Raw file', 'Intensity', 'Retention time' ]]
    print(sample_group, end="\n\n")

    #grouping data by Sequence and aggregating based on count of raw files, max for
    Intensity and cv for RT
    grouped = sample_group.groupby("Sequence")
    result = grouped.agg( {'Raw file':fun_count_rawfiles,
                          'Intensity':np.max,
                          'Retention time': func_cv}).rename(columns = {'Intensity': 'max_intensity',
                              'Retention time': 'cv_retentiontime'})

#    print("\nresult: \n",result)

```

```

#filtering peptides which are present in all Raw files
result_interm = result[result["Raw file"] == num_rawfiles ]
# print("\nresult_interm: \n",result_interm)
#selecting top 5 peptides with highest intensities and choosing peptide with lowest cv
num_peptides = result_interm.shape[0]

if num_peptides > 0:
    if num_peptides > 5:
        filter_top_n = 5
    else:
        filter_top_n = num_peptides

    print('\n filter_top_n: ',filter_top_n)
    result_interm_topn = result_interm.nlargest( filter_top_n,['max_intensity'])
    print('\nresult_interm_topn: \n', result_interm_topn)
    result_selected = result_interm_topn['cv_retentiontime'].idxmin()

    #if cv_retention time is NaN, choosing peptide with highest intensity
    result_selected = np.where( result_selected is np.NaN,
result_interm['max_intensity'].idxmax(), result_selected )
    print("\n result_selected: ",result_selected)
    selected_peptides.append(result_selected.item())

# In[13]:
print("number of selected_peptides: ",len(selected_peptides))
print("selected_peptides: ",selected_peptides)

# ## Verifying the presence of selected peptides across all Raw files

# In[14]:
df_evidence_sample = df_evidence_q2[
df_evidence_q2['Sequence'].isin(selected_peptides) ]
df_evidence_sample.groupby(['Sequence', 'Raw file']).count()['id'].unstack()

# #### Above table shows the distribution of peptides across all Raw files

# ## Info needed for each peptide:
#     Peptide sequence - To be extracted from evidence file
#     m/z           - Evidence file
#     mass error (ppm) - Evidence file

```

```

#      Retention length - Evidence file
#      Intensity      - Evidence file
#      Resolution     - Evidence file
#      1/k0          - Evidence file
#      CCS            - Evidence file

# In[15]:
evidence_columns = ['id', 'Sequence', 'Raw file', 'm/z', 'Mass error [ppm]', 'Intensity',
                    'Resolution', '1/K0', 'CCS', 'Retention time' ]
df_evidence_mod = df_evidence_sample[evidence_columns]
df_evidence_mod.head()

# ## Calculating CV and mean for Retention time, Intensity and other columns
# In[16]:
final_evidence_columns = ['Sequence', 'm/z', 'Mass error [ppm]', 'Intensity', 'Resolution',
                           '1/K0', 'CCS',
                           'Retention time']

df_evidence_mod_final =
df_evidence_mod[final_evidence_columns].groupby('Sequence').agg(
    mean_retention_time = ('Retention time', 'mean'), cv_retention_time = ('Retention
time', func_cv),
    mean_intensity = ('Intensity', 'mean'), cv_intensity = ('Intensity', func_cv),
    mean_mz = ('m/z', 'mean'),
    mean_masserror_ppm = ('Mass error [ppm]', 'mean'),
    mean_resolution = ('Resolution', 'mean'),
    mean_1_k0 = ('1/K0', 'mean'),
    mean_ccs = ('CCS', 'mean')
)
df_evidence_mod_final

# # Extracting Fragment ion info from msms data
# In[17]:
df_msms.describe()

# ### Filtering msms file with selected peptides list and charge = 2
# In[18]:
msms_columns = ['id', 'Sequence', 'Raw file', 'Matches', 'Intensities', 'Masses', 'Evidence
ID']
df_msms_filtered = df_msms.query("Sequence in @selected_peptides and Charge ==
2")[msms_columns]

```

```

df_msms_filtered.head()

# In[19]:
df_msms_filtered.describe()

# In[20]:
def fun_matches_pattern(element_matches):
    list_matches = element_matches.split(sep = ';')

    #regex pattern to detect "y/b<1 or 2 digit numbers>"
    regex = re.compile('^[y,b][0-9]{1,2}$')
    pattern_list_matches = []

    for word in list_matches:
        pattern = regex.findall(word)
        if len(pattern) == 0:
            pattern = [0]
        pattern_list_matches.append(pattern[0])

    return pattern_list_matches

# In[21]:
df_msms_filtered.loc[:, 'pattern_matches'] =
df_msms_filtered['Matches'].map(fun_matches_pattern)
df_msms_filtered.head()

# In[22]:
# function to filter masses and intensities for ions with filtered pattern
def func_filter_cols(df_input ):
    list_intensities = df_input[ 'Intensities'].split(sep = ';')
    list_masses = df_input[ 'Masses'].split(sep = ';')
    sel_list_intensities = []
    sel_list_matches = []
    sel_list_masses = []
    for index, item in enumerate(df_input['pattern_matches']):
        if item != 0:
            sel_list_matches.append(item)
            sel_list_intensities.append(list_intensities[index])
            sel_list_masses.append(list_masses[index])

    return (sel_list_matches,sel_list_intensities,sel_list_masses)

```

```

df_msms_filtered['filtered_cols'] = df_msms_filtered.apply(func_filter_cols, axis = 1)
df_msms_filtered.head()

# ## Processing filtered msms file to extract top 5 ions with highest intensity for each
sequence

# In[24]:
#function to transform matches, intensities and masses from single column to
#multiple columns
def fun_create_matches_df(input):
    print( 'input:\n',input)
    df_input = pd.DataFrame({ 'matches': input[0],
                              'intensities': input[1],
                              'masses': input[2]})
    return df_input

#create final msms file
df_msms_final = pd.DataFrame(index = df_evidence_mod_final.index)

#create dataframe with separate data for each fragment ion
df_frag_ions = pd.DataFrame(columns = ['Sequence', 'matches', 'intensities', 'masses'])

#grouping the msms file to extract fragment ion info
grouped_sequence = df_msms_filtered.groupby('Sequence')

for index, group in grouped_sequence:
    print("-" * 27)
    print('peptide sequence : ',index)
    print("-" * 27)
    group_mod = group[['Sequence', 'filtered_cols']].reset_index(drop = True)
    # print('group_mod: \n',group_mod )
    #creating dataframe to extract all fragment ions data for a peptide
    column_names = ['matches', 'intensities', 'masses']
    df_required = pd.DataFrame( columns = column_names)

    #concatenate all fragment ions info for a peptide
    for element in group_mod['filtered_cols']:
        # print('element: \n', element)
        df_input = pd.DataFrame({ 'matches': element[0],
                                  'intensities': element[1],
                                  'masses': element[2]})
        # print(df_input)

```

```

df_required = pd.concat([df_required, df_input])

df_required.reset_index(drop = True, inplace = True)
print('\nfragment ions list: \n',(df_required))

#filtering fragment ions with max intensities and creating new dataframe
df_required['intensities'] = pd.to_numeric(df_required['intensities'])
df_required['masses'] = pd.to_numeric(df_required['masses'])
df_required_with_max =
df_required.groupby('matches')['intensities'].transform(np.max)
max_intensity_idx = df_required_with_max == df_required['intensities']
df_max_intensity = df_required[max_intensity_idx]
# print('grouped df_required: \n', df_max_intensity )
df_test = df_max_intensity.drop_duplicates('matches').sort_values('intensities',
ascending=False).iloc[:5]
df_test['Sequence'] = index
print('\nchosen fragment ions: \n', df_test)
#appending fragment ion data for all peptides
df_frag_ions = df_frag_ions.append( df_test)
df_msms_final.loc[index, 'fragment_ions'] = str(df_test.to_dict(orient = 'list'))

# In[25]:
df_frag_ions.head(10)

# # Combining evidence file data with msms file

# In[26]:
df_evid_fragions_join = df_evidence_mod_final.merge(df_frag_ions, left_on =
'Sequence', right_on='Sequence' )
df_evid_fragions_join

# ## Extracting Retention length (FWHM) and Ion mobility index length (FWHM) from
allpeptides file.

# In[27]:
df_allpeptides.head()

# In[28]:
df_allpeptides.columns

```

```

# In[29]:
msms_columns = ['id', 'Sequence', 'Raw file', 'Scan number']
df_msms_test = df_msms.query("Sequence in @selected_peptides and Charge ==
2")[msms_columns]
df_msms_test.head()

# In[30]:
allpept_columns = ['Sequence', 'Raw file', 'Scan number', 'Retention length (FWHM)', 'Ion
mobility index length (FWHM)']
df_allpept_merged = df_msms_test.merge(df_allpeptides, how = 'left', left_on= ['Scan
number', 'Raw file'],
                                     right_on = ['MS/MS scan number', 'Raw file'])[allpept_columns]
df_allpept_merged.head()

# In[31]:
df_allpept_final = df_allpept_merged.groupby('Sequence').mean()[['Retention length
(FWHM)', 'Ion mobility index length (FWHM)']]
df_allpept_final

# In[32]:
df_benchmark_peptides = df_evid_fragions_join.merge(df_allpept_final, left_on=
'Sequence', right_index = True)
df_benchmark_peptides

# In[33]:
df_benchmark_peptides.head(10)

# # Creating new directory and saving final benchmarks peptides excel file

# In[34]:
# detect the current working directory and print it
current_path = os.getcwd()
print ("The current working directory is %s" % current_path)

# define the name of the directory to be created
new_path = "LCMS_data"
try:
    os.mkdir(new_path)
except OSError as error:

```



```
    print ("Creation of the directory {0} failed because of error {1}".format( new_path,
error))
else:
    print ("Successfully created the directory %s " % new_path)
finally:
    os.chdir( os.path.join(current_path, new_path))
```

```
# In[35]:
#write to excel dataset
filepath = 'Benchmark_peptides.xlsx'
with ExcelWriter(filepath) as writer:
    df_benchmark_peptides.to_excel(writer, sheet_name = ( 'df_benchmark_peptides' )
    writer.save()
```

Histograms of CV for Retention time and Intensity

```
# In[36]:
plt.figure(figsize= (20,10))
histogram_cv_RT =
sns.distplot(df_benchmark_peptides['cv_retention_time']).set_title('Histogram of CV for
Retention time')
histogram_cv_RT.get_figure().savefig('histogram_cv_RT.png')
histogram_cv_RT
```

```
# In[37]:
plt.figure(figsize= (20,10))

histogram_cv_intensity =
sns.distplot(df_benchmark_peptides['cv_intensity']).set_title('Histogram of CV for
Intensity')
histogram_cv_intensity.get_figure().savefig('histogram_cv_intensity.png')
histogram_cv_intensity
```

8.2. Part 2

```
#!/usr/bin/env python
# coding: utf-8
# # Code to extract data from test set spectrometer data based on benchmarked peptides

# In[1]:
import pandas as pd
import numpy as np
import os
from pandas import ExcelWriter
import seaborn as sns
import matplotlib.pyplot as plt
import re
import warnings

# # Please enter input file path in test_directory variable

# In[2]:
####MQoutput path
test_directory = r'/Volumes/Transcend/OneDrive/UoH DS/Summer 2020/Spectrometry lab
internship/new_MQoutput/TestDataSet-2'

# In[3]:
print ("The current working directory is %s" % os.getcwd())

####benchmark peptides path
new_path = r"/LCMS_data"
current_dir = str(os.getcwd())
if current_dir.endswith(new_path):
    working_dir = current_dir
else:
    working_dir = current_dir + new_path

#changing working directory to benchmark peptides path
os.chdir(working_dir)
print ("The new working directory is %s" % os.getcwd())

# # Reading all files in the test/txt directory
# In[4]:
```

```

for root, dirs, files in os.walk(test_directory, topdown=False):
    file_names = []
    print('root: ', root)
    print('filenames--')
    for name in files:
        if name[0] is not '.':
            print(os.path.join(root, name))
            file_names.append(os.path.join(root, name))

```

```

# In[5]:
#Reading filenames of evidence and msms files
for name in file_names:
    if name.endswith('evidence.txt'):
        file_evidence = name
    elif name.endswith('msms.txt'):
        file_msms = name
    elif name.endswith('allPeptides.txt'):
        file_allpeptides = name

```

```
warnings.filterwarnings("ignore")
```

```

# In[6]:
print(file_evidence)
df_evidence = pd.read_table(file_evidence)
df_evidence.columns = (df_evidence.columns).str.lower()
df_evidence.head()

```

```

# In[7]:
print(file_msms)
df_msms = pd.read_table(file_msms)
df_msms.columns = (df_msms.columns).str.lower()
df_msms.head()

```

```

# In[8]:
print(file_allpeptides)
df_allpeptides = pd.read_table(file_allpeptides)
df_allpeptides.columns = (df_allpeptides.columns).str.lower()
df_allpeptides.head()

```

```

# # Read extracted benchmark peptides file
# In[9]:
file_benchmark_peptides = 'Benchmark_peptides.xlsx'

```

```

df_extracted_peptides = pd.read_excel(file_benchmark_peptides, index_col=0)
df_extracted_peptides.columns = (df_extracted_peptides.columns).str.lower()

df_extracted_peptides.head()

# # Process peptides data for visualization

# In[10]:
#showing list of benchmark peptide sequences
selected_peptides = df_extracted_peptides['sequence'].unique()
selected_peptides

# ## Filter evidence and allpeptides

# In[11]:
evidence_columns = [ 'sequence', 'raw file', 'm/z', 'mass error [ppm]', 'intensity',
                    'resolution', '1/k0', 'ccs', 'retention time' ]
df_evidence_viz = df_evidence.query("sequence in @selected_peptides")[evidence_columns]
df_evidence_viz

# In[12]:
msms_columns = ['id', 'sequence', 'raw file', 'scan number']
df_msms_test = df_msms.query("sequence in @selected_peptides")[msms_columns]
df_msms_test.head()

# In[13]:
allpept_columns = ['sequence', 'raw file', 'scan number', 'retention length (fwhm)', 'ion mobility
index length (fwhm)']
df_allpept_merged = df_msms_test.merge( df_allpeptides, how = 'left', left_on= ['scan
number', 'raw file'],
                                     right_on = ['ms/ms scan number', 'raw file'])[allpept_columns]
df_allpept_merged.head()

# ## Reframing evidence file to include Rawfile name and date of processing

# In[14]:
#extracting date from Raw file name
def fun_date_extract(rawfile_names):
    #regex pattern to detect last numbers after underscore

```

```

regex = re.compile('_[0-9]+')
date_extract = regex.findall( rawfile_names)

return date_extract[-1]

# In[15]:
df_evidence_viz['date_extracted'] = df_evidence_viz['raw file'].map( fun_date_extract)
df_evidence_viz

# In[16]:
df_evidence_aggviz = df_evidence_viz.groupby( ['sequence', 'date_extracted']).mean()
df_evidence_aggviz

# ## Extracting FWHM from allpeptides file and separating date of extraction

# In[17]:
df_allpept_merged['date_extracted'] = df_allpept_merged['raw file'].map( fun_date_extract)
df_allpept_merged

# In[18]:
df_allpept_aggviz = df_allpept_merged.groupby( ['sequence', 'date_extracted']).mean()[
['retention length (fwhm)', 'ion mobility index length (fwhm)']]
df_allpept_aggviz

# ## Merging evidence and all peptides files

# In[19]:
df_evid_allpept_viz = df_evidence_aggviz.merge(df_allpept_aggviz, left_index=True,
right_index=True)
df_evid_allpept_viz

# In[20]:
df_peptides_finalviz = df_evid_allpept_viz.reset_index().copy()
df_peptides_finalviz

# In[21]:
df_peptides_finalviz.head(10)

```

```

# # Process fragment ion data for visualization

# ## Filter msms file and extract fragment ions

# In[22]:
msms_columns = ['id','sequence', 'raw file', 'matches', 'intensities', 'masses']
df_msms_filtered = df_msms.query("sequence in @selected_peptides")[msms_columns]
df_msms_filtered

# In[23]:
df_msms_filtered['date_extracted'] = df_msms_filtered['raw file'].map( fun_date_extract)
df_msms_filtered

# In[24]:
#grouping the msms file to extract fragment ion info
grouped_msms = df_msms_filtered.groupby(['sequence', 'date_extracted'])

#creating dataframe to extract all fragment ions data for a peptide
column_names = ['sequence','date_extracted','matches','intensities', 'masses']
df_fragment_ions = pd.DataFrame( columns = column_names)

#Read each group and extract fragment ions parameters and create a dataframe
for index, group in grouped_msms:
    print("-" * 27)
    print('peptide sequence: ',index)
    print("-" * 27)
    group_mod = group[['sequence', 'date_extracted','matches', 'intensities',
'masses']].reset_index(drop = True)

    #extract fragment ion parameters into list variables
    list_matches = group_mod[ 'matches'].str.split(';').tolist()[0]
    list_intensities = group_mod[ 'intensities'].str.split( ';').tolist()[0]
    list_masses = group_mod[ 'masses'].str.split(';').tolist()[0]

    req_sequence = index[0]
    benchmark_ions = list(df_extracted_peptides.query('sequence ==
@req_sequence')['matches'])
    print('benchmark_ions: \n',benchmark_ions )

    #concatenate all fragment ions info for a peptide

```

```

for ion in list_matches:
#     print('ion: \n', ion)
    if ion in benchmark_ions:
        req_index = list_matches.index(ion)
        df_input = pd.DataFrame({'sequence':index[0],
                                'date_extracted' : index[1],
                                'matches': [list_matches[req_index]],
                                'intensities': [list_intensities[req_index]],
                                'masses': [list_masses[req_index]]
                                })

        df_fragment_ions = pd.concat([df_fragment_ions, df_input])

    df_fragment_ions.reset_index(drop = True, inplace = True)
#     print('df_fragment_ions: \n', df_fragment_ions)

#convert columns to float
df_fragment_ions['intensities'] = pd.to_numeric(df_fragment_ions['intensities'])
df_fragment_ions['masses'] = pd.to_numeric(df_fragment_ions['masses'])

df_frag_ions_finalviz = df_fragment_ions

# In[25]:
df_frag_ions_finalviz.head(10)

# # Update existing peptides and ions sheets with new data

# In[26]:
current_dir = os.getcwd()
required_file = 'Peptides_viz.xlsx'

#Read all filenames in current directory and create dataframes for old data
for *rest, files in os.walk(current_dir, topdown=False):
    print('All files in current directory: ', files)
    #Read existing peptides file if it exists
    if required_file in files:
        print('FOUND REQUIRED FILE: ', required_file)
        df_peptides_finalviz_old = pd.read_excel(required_file, sheet_name= 'peptides_viz_data',
index_col=0)
        df_frag_ions_finalviz_old = pd.read_excel(required_file, sheet_name=
'fragment_ions_vizdata', index_col=0)
    else:
        print('NOT FOUND REQUIRED FILE: ', required_file)
        df_peptides_finalviz_old = pd.DataFrame( columns= df_peptides_finalviz.columns)

```

```

df_frag_ions_finalviz_old = pd.DataFrame( columns= df_frag_ions_finalviz.columns)

print('\nOld peptides_finalviz data: \n', df_peptides_finalviz_old.head())
print('\nOld frag_ions_finalviz data: \n', df_frag_ions_finalviz_old.head())

# In[27]:
#create new dataframes by combining old data with new data
df_peptides_finalviz_new = pd.concat( [df_peptides_finalviz_old, df_peptides_finalviz])
df_frag_ions_finalviz_new = pd.concat( [df_frag_ions_finalviz_old, df_frag_ions_finalviz])

print('\n updated peptides_finalviz : \n', df_peptides_finalviz_new.head())
print('\n updated frag_ions_finalviz : \n', df_frag_ions_finalviz_new.head())

# In[28]:
#write new files into excel file
filepath = 'Peptides_viz.xlsx'
with ExcelWriter(filepath) as writer:
    df_peptides_finalviz_new.to_excel(writer, sheet_name = ( 'peptides_viz_data' ) )
    df_frag_ions_finalviz_new.to_excel(writer, sheet_name = ( 'fragment_ions_vizdata' ) )
    writer.save()

```

8.3. Part 3

```

#!/usr/bin/env python
# coding: utf-8

# # Code to plot peptides and fragment ions visualization

# In[1]:
import pandas as pd
import numpy as np
import os
from pandas import ExcelWriter
import cufflinks as cf
import plotly.offline
import warnings

# # Read Peptides data file
# In[2]:
print ("The current working directory is %s" % os.getcwd())

####Processed peptides data path
working_dir = r'/Users/kishore/LCMS_data'
#changing working directory to above path
os.chdir(working_dir)

```



```

#Reading the file
peptides_finalviz = 'Peptides_viz.xlsx'
df_peptides_finalviz = pd.read_excel(peptides_finalviz, sheet_name='peptides_viz_data', index_col= 0)
df_frag_ions_finalviz = pd.read_excel(peptides_finalviz, sheet_name='fragment_ions_vizdata',
index_col= 0)
print('df_extracted_peptides: \n', df_peptides_finalviz.head())
print('df_frag_ions_finalviz: \n', df_frag_ions_finalviz.head())

# # Data visualizations
# ## Visualization of date vs parameters using plotly library

# In[3]:
cf.go_offline()
cf.set_config_file(offline=True, world_readable=True)
plotly.offline.init_notebook_mode()

# In[4]:
warnings.filterwarnings("ignore")
for col in df_peptides_finalviz.columns[2:]:
    df_peptides_finalviz.iplot(kind = 'scatter', mode = 'lines+markers',
        x = 'date_extracted', y = str(col), categories = 'sequence',
        title = 'date vs '+str(col), xTitle = 'date', yTitle = str(col))

# # Extracting Fragment ion info from msms data
# In[5]:
df_frag_ions_finalviz.head()

# In[6]:
warnings.filterwarnings("ignore")

#showing list of benchmark peptide sequences
selected_peptides = df_frag_ions_finalviz['sequence'].unique()

for peptide in selected_peptides:
    df_peptide = df_frag_ions_finalviz.query('sequence == @peptide')[['sequence', 'date_extracted',
'matches', 'intensities']]
    print('Cross table for above peptide: \n', df_peptide)
    if df_peptide.empty:
        print('\n#####No fragment ion data for this peptide#####\n')
    else:
        df_peptide.iplot(kind = 'scatter', mode = 'lines+markers',
            x = 'date_extracted', y = 'intensities', categories = 'matches',
            title = 'date vs fragment ion intensities for '+peptide, xTitle = 'date', yTitle = 'intensities')

```